

MANUAL CONTROLLER RD-SC11



Product Specification

MANUAL CONTROLLER RD-SC11

Product Characteristics

- This manual controller is oriented to the elderly care bed market. It is specially designed according to the physical characteristics of the middle-aged and elderly, and humanized design is carried out from appearance, color matching to function; Bright color matching, comfortable grip, and absorbing the product advantages of Runda's many other medical systems in terms of functions, refining and upgrading to create a new type of hand controller.



Product Characteristics

- Protection grade: IPX6
- Number of keys: Up to 16 keys can be designed
- Drive up to 5 actuators simultaneously
- Color: grayish white
- This manual controller is designed for the elderly. The product uses silicone buttons, which are comfortable to touch
- Removable manual controller cable
Controller kit: Runda all conventional controllers
- Both hands can operate easily
- Panel color and style can be customized
- The shell feels comfortable, with soft lines and beautiful appearance

USE

- Operating temperature: 5-40 °C
- Storage temperature: - 5-45 °C
- Supporting controller: optional Runda conventional controller
- Relative humidity: 20% - 80% at 30 °C, non condensing
- Atmospheric pressure: 860-1060hPa
- Altitude: up to 2000m
- Safety certification: comply with ROHS specification

Authentication


- EN60601-1/EN60601-1-2 & ROHS2.0

Naming Rules

Model	RD-SC11	
Number of key rows (Matched with controller)		1=1row 2=2rows 3=3rows 4=4rows 5=5rows 6=6rows 7=7rows 8=8rows
Expand functions		A=message B=backlight C=battery D=light strip F=reset J=memory K=startup L=linkage R=heating S=locked T=emergency stop Y=voice Z=vibration 0=no function
Colour		G=Grey B=Black
Wire plug		Z=Straight line W=Bending
Indicator light		1=With indicator 2=Without indicator
Spring wire length*D*		1=natural length: 1m (stretching length: 3m) D=customized for special requirements
Linkage mode		0=All push rods work independently L=Multiple push rods work simultaneously (icon shall prevail)


RD-SC11 - 8 - - - G - Z - 1 - D - 0

示例




RD-SC11 standard panel (facial mask icon can be customized)

When ordering, please check whether the manual function corresponds to the symbol on the panel.



Product Dimension Drawing



Wire plug

Example

In-line



*Note:

The number of matching push rods is different, the model of matching controller is different, and the plug code will be different. Before placing an order, please contact the salesperson in the corresponding area of Runda for confirmation.

Conditions of Use

Users must determine the applicability of Runda products under specific conditions. Runda is committed to providing accurate and up-to-date product information. However, due to the continuous improvement of Runda's products, Runda's products may be frequently changed without prior notice. Therefore, Runda cannot guarantee the correctness and authenticity of the product information. Although Runda will try its best to meet the order requirements, for the above reasons, Runda cannot guarantee the delivery of any specific product. Therefore, Runda reserves the right to stop selling any products listed in the website or product catalog or other written materials of Runda. All sales activities must be carried out in accordance with the requirements of Jiechang Standard Terms of Sales and Delivery. For copies of relevant documents, please contact Runda.